



# Effects of Sterilization and 1-Year Real-Time Aging on Material Tensile Strength (MD) for 1059B



**Transition Protocol material performance is equivalent to, or better than, current Tyvek®.**

# Effects of Sterilization and 1-Year Real-Time Aging on Material Tensile Strength (CD) for 1059B




- ▲ Transition Protocol Material 1-Year Real-Time Aging
- ▲ Control 1-Year Real-Time Aging
- Transition Protocol Material 0-Year
- Control 0-Year

ASTM D5034  
 Control = DuPont™ Tyvek® 1059B  
 Center point = 0 lb<sub>f</sub>/4 in.  
 Outer point = 150 lb<sub>f</sub>/4 in.  
 CD = Cross Direction

**Transition Protocol material performance is equivalent to, or better than, current Tyvek®.**

# Effects of Sterilization and 1-Year Real-Time Aging on Material Elongation (MD) for 1059B




**Transition Protocol material performance is equivalent to, or better than, current Tyvek®.**

# Effects of Sterilization and 1-Year Real-Time Aging on Material Elongation (CD) for 1059B



Tyvek.




- ▲ Transition Protocol Material 1-Year Real-Time Aging
- ▲ Control 1-Year Real-Time Aging
- Transition Protocol Material 0-Year
- Control 0-Year

ASTM D5034  
 Control = DuPont™ Tyvek® 1059B  
 Center point = 0%  
 Outer point = 30%  
 CD = Cross Direction

**Transition Protocol material performance is equivalent to, or better than, current Tyvek<sup>®</sup>.**

# Effects of Sterilization and 1-Year Real-Time Aging on Material Puncture Strength for 1059B



- ▲ Transition Protocol Material 1-Year Real-Time Aging
- ▲ Control 1-Year Real-Time Aging
- Transition Protocol Material 0-Year
- Control 0-Year

ASTM F1342  
 Control = DuPont™ Tyvek® 1059B  
 Center point = 0 lb<sub>f</sub>  
 Outer point = 3.5 lb<sub>f</sub>

**Transition Protocol material performance is equivalent to, or better than, current Tyvek®.**